




Produced by JCM Seating Solutions Ltd  
Tel: 01733 405830 Email: enquiries@jcmseating.co.uk Web: www.jcmseating.co.uk

# Proxim<sup>®</sup>at

## User Manual



A simple, child friendly method of measuring hip range movement in children.

# CONTENTS

Background & Reliability of the Proximat.....	1
What is the Proximat? .....	1
The Components of the Proximat .....	1
How is the Proximat Used? .....	2
Note - Before Use:.....	2
Additional Measuring Facility .....	2
How is the Hip Joint Aligned? .....	3
To Check the Alignment / Starting Position of the Hip.....	3
To Measure Abduction.....	5
To Measure Adduction.....	6
To Measure Internal Rotation .....	7
To Measure External Rotation .....	8
Example of Use for Younger Children.....	9
References .....	10

## Background of the Proximat

The Proximat was developed by physiotherapists in clinical practice. It is a simple child friendly tool for measuring the hip range of movement in children. It is aimed at overcoming pragmatic difficulties with goniometry.

A service development introduced new postural management programmes aimed at preventing deformity pain and functional difficulties<sup>1-6</sup>. As a part of monitoring this programme physiotherapists were regularly measuring hip range of movement. However this was difficult with conventional goniometry and the results appeared inaccurate. Subsequent searching of the literature revealed goniometry for children with cerebral palsy has limited reliability<sup>7-12</sup>.

## Reliability of the Proximat

Following successful trial use within the service a study was carried out to assess the reliability of the proximat. Passive hip abduction, adduction, external and internal rotation were measured using the Proximat on 26 children with cerebral palsy. Testing was performed by two physiotherapists to assess inter-rater reliability and then repeated the following day to assess test-retest reliability. The Proximat was quick and easy to use and acceptable to the children. High reliability was found for all movements (ICCs = 0.83-0.93)<sup>13</sup>.

## The components of the Proximat are:

- Flexible mat
- Dry wipe pen
- Tape measure
- Recording sheet
- Poster
- CD containing all of the Proximat literature and record sheet.
- Protective case

## How is the Proximat Used?

The mat is placed on the treatment bed or floor. The child is placed on the mat (hoisted if necessary) so that the hip joint lies over the centre point + of the mat and the leg lies along the vertical 0° axis with pelvis level.

Before measuring, any excess muscle tone should be inhibited as much as possible and the child's head positioned in mid line.

The pelvis is stabilised by giving comfortable downward pressure over the iliac crests while the leg is being moved.

## Note - Before Use:

The Proximat is intended for use by a qualified clinician. Assistance may be given by anyone acting under their direction.

## Additional Measuring Facility

Any bony protrusions or the child's entire height can be measured by marking points on the mat with dry wipe pen and measuring with the tape measure.

This pen can also be used to temporarily record the measurements on the mat before transferring them to case notes.

## Proximat Safety Notice

Please note that due to the design of the Proximat it has a very smooth surface. Users are to be aware that the surface could be slippery, especially if walked on in socks. It is the users responsibility to ensure that all parties are aware of this and take due care and attention to avoid accidents.

## How is the hip joint aligned?

The anterior aspect of the hip joint lies approximately 1.5 cms below the mid point of an imaginary line between the Anterior Superior Iliac Spine and the Pubic Symphysis.




## To check the alignment / starting position of the hip

(Refer to diagrams opposite, diagrams show example of right leg positioning)


When the leg is aligned with the vertical axis 0° on the mat, the accuracy of the hip location can be confirmed by moving the leg out to the side and checking that it remains aligned to the centre point + of the Proximat.

If the hip comes out of alignment, it should be repositioned on the centre point +.



When the hip is not aligned the leg crosses the radial lines when moved:



Vertical axis 0° Line




If this happens, re-align the hip to the centre point:



Centre (+) Point

Example of hip joint aligned to the centre point:



## To Measure Abduction

- Stabilise the pelvis.
- Move the leg out to the side, away from the body.
- Stop when end of range is reached.
- Read and note the angle achieved from the mat.



## To Measure Adduction

- Stabilise the pelvis.
- Raise the opposite leg.
- Move the leg to be measured under the raised leg.
- Stop when end of range is reached.
- Read and note the angle achieved from the mat.



## To Measure Internal Rotation

- Stabilise the pelvis.
- Bend the leg to be measured to 90° at the hip and knee.
- Keep the long axis of the femur directly over the hip joint (and centre point + of the Proximat)
- Rotate the leg so the foot moves away from the other leg.
- Stop when end of range is reached.

### Take the reading from the Proximat

The degree of internal rotation is read from the bent leg's alignment with the markings on the mat as viewed from directly above.



## To Measure External Rotation

- Stabilise the pelvis.
- Bend the leg to be measured to 90° at the hip and knee.
- Keep the long axis of the femur directly over the hip joint (and centre point + of the Proximat)
- Rotate the leg so that the foot moves over the other leg
- Stop when end of range is reached.

### Take the reading from the Proximat

The degree of external rotation is read from the bent leg's alignment with the markings on the mat as viewed from directly above.



## Example of Use for Younger Children

The example below shows measurement of right hip abduction in a young child. The pelvis has been stabilised by the clinician doing the measurement.



## References

1. Gericke T. (2006) Postural Management for children with Cerebral Palsy: consensus statement. *Dev Med Child Neurol*, Vol 48, no 4: 244
2. Scrutton D, Baird G. (1997) Surveillance measures of the hips of children with bilateral cerebral palsy. *Arch Dis Child*, vol. 76, no.4: 381-384.
3. Dobson F, Boyd RN, Parrott J, Nattrass J, Graham HK. (2002). Hip Surveillance in Children with cerebral palsy: impact on the management of spastic hip disease. *J Bone Joint Surg*, vol.84B, no 5: 720-726.
4. Hagglund G, Andersson S, Duppe H, Lauge-Pederson H, Nordmark E, Westbom L. (2005) Prevention of dislocation of the hip in children with cerebral palsy. The first ten years of a population based prevention programme. *J Bone Joint Surg Br*, vol.87,no.1:95-101.
5. Scrutton D. (1989) The early management of hips in cerebral palsy. *Dev Med Child Neurol*, vol.31(1):108-116.
6. McKearnan KA, Kiechhefer GM, Engel JM, Jensen MP, Labyak S.(2004) Pain in children with cerebral palsy: a review. *J Neurosci Nurs*,vol.36 no.5:252-259.
7. Cadenhead SL, McEwen IR, Thompson DM. (2002) Effect of passive range of motion exercises on lower- extremity goniometric measurements of adults with cerebral palsy: a single subject design. *Phys Ther*, vol 82,no7,:658-659.
8. Cadenhead SL, McEwen IR, Thompson DM. (2002) Effect of passive range of motion exercises on lower- extremity goniometric measurements of adults with cerebral palsy: a single subject design. *Phys Ther*, vol 82,no7,:658-659.
9. Ashton BB, Pickles B, Roll JW. (1978) Reliability of goniometric measurements of hip motion in spastic cerebral palsy. *Dev Med Child Neurol*, vol 20 no.1 :87-94.
10. Stuberg WA, Fuchs RH, Miedaner JA. (1988). Reliability of goniometric measurements of children with cerebral palsy. *Dev Med Child Neurol*, vol 30, no. 5, :657-666.
11. McDowell BC, Hewitt V, Nurse A, Weston T, Baker R. (2000) The variability of goniometric measurements in ambulatory children with spastic cerebral palsy. *Gait Posture* vol.11(2):86-91.
12. Fosang AL, Galea MP, McCoy AT, Reddihough DS, Stofy I. (2003) Measures of muscle and joint performance in the lower limb of children with cerebral palsy. *Developmental Medicine and Child Neurology*. Vol 45, no 10, 664-670
13. Pott P, Selley A, Tyson S. (2007). The reliability and measurement error of the Proximat: a new tool for measuring hip range of movement in children with cerebral palsy. *APCP conference 2007 Free paper presentation*.