



# BeSafe Car Seat



## BeSafe iZi Up

Suitable from around 5 -12 years / 15-36 kg

### Standard features for BeSafe iZi Up

- Crash tested with all accessories
- Special headrest, offering variable support
- Extra thigh support
- Side wings, adjustable in breadth and height
- Pommel
- Soft seat unit
- Soft back padding
- Adjustable seat depth and back height
- Adjustable back angle
- Ready for fitting turntable
- Ready for fitting vest
- Ready for crotch strap

BeSafe iZi UP is now available for children and young people up to around 135 cm with stability problems, who need more support and comfort. The new solution has been tested with adjustable side wings, pommel, adjustable headrest, vest, crotch strap, turntable and extra thigh support. All tested in accordance with ECE4404.



If you have problems with the child sliding out of the seat unit, you can choose the crotch strap solution.

The headrest provides optimum support. A great solution for the users and their safety.

The pommel is easy to remove.

iZi Up can be supplied with a vest retainer, to hold the child in place during the car ride. iZi Up is supplied in dark grey.

The turntable allows you to turn the seat out towards the car door and put the user in the seat, before turning it in again. You can hear a click when the turntable has locked in the direction of travel. It may be necessary for larger users to have the back supported while the seat is turned.

Both the angle and height of the back can be adjusted. The back angle can be adjusted by around 15 degrees by the lever at the side of the seat. The back returns automatically to the upright position when the child is removed from the seat.



# www.jcmseating.co.uk



JCM Seating Solutions  
15-18 Maxwell Road  
Woodston Industrial Estate  
Peterborough  
Cambridgeshire  
PE2 7HU  
United Kingdom

**T:** 01733 405830  
**F:** 01733 405838  
**E:** enquiries@jcmseating.co.uk  
**W:** www.jcmseating.co.uk