

# TX Size 1 Specifications



Make & Model	TX-30 S1
Width (Interface with Frame)	330
Overall Width (Castor center)	490
Wheelbase	372 - 447
O/L Length	588 - 663
Minimum Turning Radius	500
Interface Height	440
Push Handle Height	560 - 1110
Tilt Angle	+5/-30
Front Wheels	150 (6")
Rear Wheels	300 (12")/ 400 (16"0)/ 600 (24")
Chassis Weight	15Kg Provisional
Max. User Weight	75Kg
Tested	ISO 7176 PART 19
Type	Att. Prop
Chair Selection (Guideline)	T0 – T2 / NEP1