

StarX Fabric Options

All JCM fabrics are fire retardant and have high standards in durability and stain resistance. If you require any further assistance our national team of representatives are always on hand.



Stirata - Turquoise



Stirata - Hot Pink

Boreal (Vinyl)

A vinyl with a unique soft pearlescent finish. Its subtle iridescent sheen is faintly metallic, which is both lustrous and warm to the touch. Boreal is tough and durable, making it ideal for heavy usage.

Cleaning and Care


Regularly wipe over with a damp soapy cloth. Clean stains immediately to avoid absorption into the vinyl. Strong detergents, solvents, polish, wax or patent chemical cleaners must not be used as they could damage the material.

- Tough & durable (non machine washable)
- Suitable for heavy usage.
- Fireproof: BS 5852: 1990 Ignition source 5.



Frame Colours

You can choose to have either flat or translucent colour (depending on model of chair chosen, an upcharge may apply). Actual colours may vary slightly.



Stirata (Soft and Flexible)

A revolutionary, high performance 4-way stretch upholstery fabric which is soft to the touch and very flexible.

- Anti-bacterial/Anti-fungal
- Waterproof: Hydrostatic Head BS3424 > 1 metre.
- Breathable: Moisture vapour permeable BS3546 65%
- Martindale Abrasion: End point > 100,000 – No pilling
- Tear Strength: BS4303 Severe Contract Grade
- Tensile Strength: BS2576
- 1972 Severe Contract Grade
- Extremely Stain Resistant
- Flame Retardant: BS5852 Crib 5

Cleaning and Care


Wipe frequently with a damp cloth to maintain the appearance of the fabric. Most stains can be removed with cold/warm

water and a clean cloth. Use cold water when treating stains caused by bodily fluids. If soap based cleaners are used, all soap residues must be rinsed from the

fabric. Can be cleaned with mild hypochlorite solutions (half a cup of household bleach with 5 litres of water) or mineral spirit. Can be steam cleaned in situ.



Stirata Illusions (mottled effect)



PLEASE NOTE! Due to printing techniques and screen previews the colours of fabrics shown on this brochure may differ slightly from the actual fabrics. We recommend that you request a fabric swatch or view the fabric booklet via your sales representative prior to placing your order.