

TX Size 1 Specifications



Make & Model	TX-30 S1
Width (Interface with Frame)	330
Overall Width (Castor center)	490
Wheelbase	372 - 447
O/L Length	588 - 663
Minimum Turning Radius	500
Interface Height	440
Push Handle Height	560 - 1110
Tilt Angle	+5/-30
Front Wheels	150 (6")
Rear Wheels	300 (12")/ 400 (16"0/ 600 (24")
Chassis Weight	15Kg Provisional
Max. User Weight	75Kg
Tested	ISO 7176 PART 19
Type	Att. Prop
Chair Selection (Guideline)	T0 – T2 / NEP1

TX Wheelchair Chassis 2011



Part	Part Number	Price
TX Size 1 chassis without interface	TX-30	\$3,623.76
Standard Features		
Drum brakes		
Adjustable length chassis		
Tilt in space 30"		
Height adjustable push handles with angle inclination		
Front wheel/castors		
6" PU Tyre		Inc
Rear wheel (choose one option)		Inc
12" PU Tyre		Inc
16" PU Tyre		Inc
Interface Options		
Standard JCM interface	TX-INT-01	Inc
Accessories		
Anti tip	TX-AT-01	Inc
Tip lever	TX-TL-01	Inc