



Fabrics



All JCM fabrics are fire retardant and have high standards in durability and stain resistance.

Stirata (Soft and Flexible)

A revolutionary, high performance 4-way stretch upholstery fabric which is soft to the touch and very flexible.

- Anti-bacterial/Anti-fungal
- Waterproof: Hydrostatic Head BS3424 > 1 metre.
- Breathable: Moisture vapor permeable BS3546 65%
- Martindale Abrasion: End point > 100,000 – No pilling
- Tear Strength: BS4303 Severe Contract Grade
- Tensile Strength: BS2576 1972 Severe Contract Grade
- Extremely Stain Resistant
- Flame Retardant: BS5852 Crib 5



Purple



Black




Hot Pink



Royal



Red



Beige

Please note that due to printing techniques and screen previews the colors shown may differ slightly from the actual fabrics. If you require any further assistance please call.

Cleaning and Care

Wipe frequently with a damp cloth to maintain the appearance of the fabric. Most stains can be removed with cold/warm water and a clean cloth.

Use cold water when treating stains caused by bodily fluids. If soap based cleaners are used, all soap residues must be rinsed from the fabric.

Can be cleaned with mild hypochlorite solutions (half a cup of household bleach with 5 liters of water) or mineral spirit. Can be steam cleaned in situ.

Moonlite Sleep System

All fabrics on our Moonlite sleep system can be machine washed at 60°C.

www.jcmamerica.comTM



JCM AmericaTM
521- D Uwharrie Court
Raleigh, NC 27606

Phone 919-831-9000
Fax 919-831-9100
Toll Free 888-831-9831
Email enquiries@jcmamerica.com
Website www.jcmamerica.com