




JCM Triton Seating Unit Wheelchair and Buggy Base Seamless Integration



We can now provide our specialist Triton™ seating system integrated seamlessly to a range of indoor and outdoor mobility base options.



Our Triton Seating System comes with a manufacturer's 3 Year Warranty

The Triton seating system comes with a choice of seat, back and body supports as well as a comprehensive list of accessories and a choice of pressure relieving cushions and fabric

options. Individually tailored support and growth adjustment means the Triton seating unit maximises positional functionality and user versatility.



quick & easy transfer

TRITON™ interfaced with the *Otto Bock*® Kimba Spring Buggy Base

The Kimba is an extremely manoeuvrable and lightweight buggy base. The spring system is specified to the weight of the child enabling a smooth, jolt-free ride. The base is made from high-quality aluminium with folding guards.

The Triton seating unit effortlessly separates from the buggy base, allowing the lightweight (15kg) base and interface to be folded down for transportation.



The Kimba Spring interfaced with the Triton combines high functionality with outstanding design and practicality.

Kimba buggy base with Triton 0 seating unit

TRITON™ interfaced with the Otto Bock Discovery Wheelchair Base

The Discovery is a particularly stable and versatile mobility base. It is suitable for both indoor

and outdoor use and is notable for its high load capacity (up to 110 kg).

The Otto Bock base is available with a range of wheel and handle options and the seat can tilt in space up to 35°.



Otto Bock Discovery with 24" wheels



Otto Bock Discovery interfaced with Triton 3 seating unit

Otto Bock®

TRITON™ interfaced with the Spectra Blitz power base

The Invacare Spectra Blitz is a versatile indoor / outdoor power chair for children. Its compact design provides excellent mobility and driving performance.

This unique combination offers the highest level of postural management with the option of either user or attendant controlled mobility.



INVACARE



The Spectra is an easily manoeuvrable chair, equipped with rear suspension for a comfortable ride.

www.jcmamerica.com™

JCM America™ 521- D Uwharrie Court, Raleigh, NC 27606

Phone: 919-831-9000 Fax: 919-831-9100 Toll Free: 888-831-9831 Email: enquiries@jcmamerica.com